

Stainless steel embedding molds

IVD In vitro diagnostic – medical device

CND number: W05030599

General features

Varied measures steel tiles for inclusion of histological paraffin samples.

| Code | Model | Size (mm) | Packaging |
|------------|-----------------|--------------|-----------|
| 07-BM15155 | Bio Mold | 15 x 15 x 5 | 12 pcs |
| 07-BM24245 | Bio Mold | 24 x 24x 5 | 12 pcs |
| 07-BM30245 | Bio Mold | 30 x 24 x 5 | 12 pcs |
| 07-BM37245 | Bio Mold | 37 x 24 x 5 | 12 pcs |
| 07-BM775 | Bio Mold | 7 x 7 x 5 | 12 pcs |
| 07-MBM6 | Mega Mold | 33 x 24 x 12 | 6 pcs |
| 07-7010 | Super Mega Mold | 65 x 45 x 15 | 10 pcs |



Manufacturer: **Bio-Optica Milano S.p.A.**

Functional features

Material: stainless steel.

Bio Molds are reusable products. To guarantee the cleaning, we advise to clean them with xylene or d-limonene, after use, to eliminate all paraffin remnants.

Validity date: not applicable.

For the characteristics of the product, particular precaution for handling and storage are not necessary.

The embedding moulds are packed in a cellophane bag, enough for guaranteeing integrity and cleaning until use.

