

for
IMMUNOSTAINING

LIQUID BLOCKER

SUPER PAP PEN

Pap-Pen has been renewed by Liquid Blocker (Super Pap Pen) with great improvement !



PAP PEN minimize the consumption of liquid for immunostaining by using this PAP PEN as a "fence" on slide glass.

It is very easy to use.
You only write a "circle" around tissue by this pen on slide glass.
Then, a fence are created on the slide glass to prevent overflow of liquid for immunostaining.

Two pen tip sizes,
4mm(Regular)
and 2mm (Mini).



MADE IN JAPAN

